

BALANZA ELECTRÓNICA - 1000G A 5000G X 0.01G



Balanza Electrónica

- Capacidad Max :1000g a 5000g
- Legibilidad : 0.01g
- Peso Bruto : 3 kg

SKU: N / A | **Categorías:** [0.01g/0.02g/0.05g](#) |

VARIACIONES

Imagen	SKU	Descripción	2021 Peso máximo (gramos)
	B-01-02-03-0301	1000	
	B-01-02-03-0302	2000	
	B-01-02-03-0303	3000	
	B-01-02-03-0304	4000	
	B-01-02-03-0305	5000	

DESCRIPCIÓN DEL PRODUCTO

Especificaciones Técnicas

Modelo	XU10002DG	XU20002DG	XU30002DG	XU40002DG	XU50002DG
Código de Producto	B-01-02-03-0301	B-01-02-03-0302	B-01-02-03-0303	B-01-02-03-0304	B-01-02-03-0305
Capacidad Max. (g)	1000g	2000g	3000g	4000g	5000g
Legibilidad (g)	0.01g				
Pantalla	LCD (luz blanca de fondo)				
Tamaño del platillo (mm)	Φ160mm				
Energía	110V/220V Opcional: batería recargable				
Unidades	g, ct, oz, lb, GN, %, piezas				
Dimensión del esquema (cm)	33×21×10				
Tamaño de embalaje (cm)	33×28×15				
Interfaz	RS232				
Opcional	Gancho para peso inferior , Batería recargable				
Peso Bruto (kg)	3				

INFORMACIÓN ADICIONAL

COTECNO