


REGULATOR OUTPUT
PRESSURE

PRESSURE
REGULATOR
ASSEMBLY
P/N 1144026

TANK PRESSURE



PRESSURE RELIEF
BLEED VALVE
ASSEMBLY
P/N 1155893

HIGH PRESSURE HOSE,
6 Ft - 3000 PSI
P/N1151427

NOTE:
SWAGELOK TIGHTENING
INSTRUCTIONS ARE ATTACHED.
PLEASE FOLLOW AS THEY
ARE VERY IMPORTANT.

ATTACH TO
HP4750

STERLITECH CORPORATION
22027 - 70TH AVE S. KENT, WA 98032

HP 4750 PRESSURE
REGULATOR
ASSEMBLY

DO NOT SCALE DRAWING