

**PROBETA DE VIDRIO BELLTRONIC DE 250 ML CON BASE HEXAGONAL DE VIDRIO,  
BELLPROBVID-250ML**



**SKU:** B-01-47-0106-1005 | **Categorías:** [Insumos de Laboratorio](#), [Material de vidrio de laboratorio](#), [Probeta de vidrio](#) |

## DESCRIPCIÓN DEL PRODUCTO

Codigo	Código SKU	Capacidad (ml)	Sub. División (ml)	Tolerancia + ml	Piezas/Contendor
BellProbVid-5ml	B-01-47-0106-1000 5		0.1	0.1	360
BellProbVid-10ml	B-01-47-0106-1001 10		0.2	0.2	288
BellProbVid-25ml	B-01-47-0106-1002 25		0.5	0.5	288
BellProbVid-50ml	B-01-47-0106-1003 50		1	1	144
BellProbVid-100ml	B-01-47-0106-1004 100		1	1	96
BellProbVid-250ml	B-01-47-0106-1005 250		2	2	48
BellProbVid-500ml	B-01-47-0106-1006 500		5	5	24
BellProbVid-1000ml	B-01-47-0106-1007 1000		10	10	12
CWGPV00001000	CWGPV00001000 1000		-	-	-
BellProbVid-2000ml	B-01-47-0106-1008 2000		20	20	8

## INFORMACIÓN ADICIONAL

COTECNO