

CONJUNTO DE PRESURIZACIÓN BELLTRONIC, VÁLVULA, 1/8 PULG.



SKU: B-01-16-02-01-1000 | **Categorías:** [Hornos de Celdas / Rodillos de Envejecimiento](#) |

DESCRIPCIÓN DEL PRODUCTO

[vc_row type="in_container" full_screen_row_position="middle" column_margin="default" scene_position="center" text_color="dark" text_align="left" overlay_strength="0.3" shape_divider_position="bottom" bg_image_animation="none"][vc_column column_padding="no-extra-padding" column_padding_position="all" background_color_opacity="1" background_hover_color_opacity="1" column_link_target="_self" column_shadow="none" column_border_radius="none" width="1/1" tablet_width_inherit="default" tablet_text_alignment="default" phone_text_alignment="default" overlay_strength="0.3" column_border_width="none" column_border_style="solid" bg_image_animation="none"][vc_column_text]**Conjunto de Presurización Belltronic, Válvula, 1/8 pulg. B-01-16-02-01-1000**[/vc_column_text][vc_column][vc_row type="in_container" full_screen_row_position="middle" column_margin="default" scene_position="center" text_color="dark" text_align="left" overlay_strength="0.3" shape_divider_position="bottom" bg_image_animation="none"][vc_column column_padding="no-extra-padding" column_padding_position="all" background_color_opacity="1" background_hover_color_opacity="1" column_link_target="_self" column_shadow="none" column_border_radius="none" width="1/1" tablet_width_inherit="default" tablet_text_alignment="default" phone_text_alignment="default" overlay_strength="0.3" column_border_width="none" column_border_style="solid" bg_image_animation="none"][vc_column_text]**Descripción del Conjunto de Presurización Belltronic, Válvula, 1/8 pulg. B-01-16-02-01-1000:**

La válvula de presurización está diseñada para usarse con celdas de envejecimiento presurizadas. Válvula de presión de acero inoxidable 316, certificada para 2500 psi (17.237 kPa)

Especificaciones de la celda de envejecimiento

Modelo	Material	Volumen de celda	Max. Presión de trabajo		Max. Temperatura.		Información de envío
			psig	kPa	°F	°C	
B-01-16-02-01-0200	Hastelloy u otras aleaciones disponibles	500ml o 260ml	3200	22000	850	450	10kg
B-01-16-02-01-0500	Acero inoxidable 304	500 mL	2500	17237	500	260	8kg
B-01-16-02-01-0700	Acero inoxidable 316	500mL	2500	17237	600	320	10kg
B-01-16-02-01-0600	Liner de teflón con celdas de envejecimiento de 500 ml	500mL	2500	17237	600	320	11kg
B-01-16-02-01-0800	Acero inoxidable 304	260 mL	2500	17237	500	260	5kg
B-01-16-02-01-0900	Acero inoxidable 316						

INFORMACIÓN ADICIONAL

COTECNO