

EATON LOFMEM WB LMWB-10-0.45-DOES FILTER CARTRIDGE



SKU: LMWB-10-0.45-DOES | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

DESCRIPCIÓN DEL PRODUCTO

EATON LOFMEM WB LMWB-10-0.45-DOES FILTER CARTRIDGE

PART NUMBER: LMWB-10-0.45-DOES

[Please contact us for more information and/or a quotation request.](#)

Specifications

Brand	Eaton LOFMEM WB
Model	LMWB 10 0.45 DOES
Pore Size μ	0.45

Dimensions / Weights

Length in inch	10.0
Length in cm	25.4
Quantity per case	12 pieces



[Download Datasheet](#)

INFORMACIÓN ADICIONAL

url_fab <https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-10-0.45-DOES/LMWB-10-0.45-DOES-Filter-Cartridge/index.html>

brand Lenntech

COTECNO