

EATON LOFMEM WB LMWB-20-1.2-4S FILTER CARTRIDGE



SKU: LMWB-20-1.2-4S | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

DESCRIPCIÓN DEL PRODUCTO

EATON LOFMEM WB LMWB-20-1.2-4S FILTER CARTRIDGE

PART NUMBER: LMWB-20-1.2-4S

[Please contact us for more information and/or a quotation request.](#)

Specifications

Brand	Eaton LOFMEM WB
Model	LMWB 20 1.2 4S
Pore Size μ	1.2

Dimensions / Weights

Length in inch	20.0
Length in cm	50.8
Quantity per case	12 pieces



[Download Datasheet](#)

INFORMACIÓN ADICIONAL

url_fab <https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-20-1.2-4S/LMWB-20-1.2-4S-Filter-Cartridge/index.html>

brand Lenntech

COTECNO