

EATON LOFMEM WB LMWB-30-1.2-3S FILTER CARTRIDGE



SKU: LMWB-30-1.2-3S | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

COTECNO

DESCRIPCIÓN DEL PRODUCTO

EATON LOFMEM WB LMWB-30-1.2-3S FILTER CARTRIDGE

PART NUMBER: LMWB-30-1.2-3S

[Please contact us for more information and/or a quotation request.](#)

Specifications

Brand	Eaton LOFMEM WB
Model	LMWB 30 1.2 3S
Pore Size μ	1.2

Dimensions / Weights

Length in inch	30.0
Length in cm	76.2
Quantity per case	12 pieces



[Download Datasheet](#)

INFORMACIÓN ADICIONAL

url_fab <https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-30-1.2-3S/LMWB-30-1.2-3S-Filter-Cartridge/index.html>

brand Lenntech

COTECNO