

EATON LOFMEM WB LMWB-30-1.2-3S FILTER CARTRIDGE



SKU: LMWB-30-1.2-3S | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

COTECNO

DESCRIPCIÓN DEL PRODUCTO

EATON LOFMEM WB LMWB-30-1.2-3S FILTER CARTRIDGE

PART NUMBER: LMWB-30-1.2-3S

[Please contact us for more information and/or a quotation request.](#)

Specifications

| | |
|-----------------|-----------------|
| Brand | Eaton LOFMEM WB |
| Model | LMWB 30 1.2 3S |
| Pore Size μ | 1.2 |

Dimensions / Weights

| | |
|-------------------|-----------|
| Length in inch | 30.0 |
| Length in cm | 76.2 |
| Quantity per case | 12 pieces |



[Download Datasheet](#)

INFORMACIÓN ADICIONAL

url_fab <https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-30-1.2-3S/LMWB-30-1.2-3S-Filter-Cartridge/index.html>

brand Lenntech

COTECNO