

## EATON LOFMEM WB LMWB-30-1.2-DOES FILTER CARTRIDGE



**SKU:** LMWB-30-1.2-DOES | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

## DESCRIPCIÓN DEL PRODUCTO

### EATON LOFMEM WB LMWB-30-1.2-DOES FILTER CARTRIDGE

**PART NUMBER: LMWB-30-1.2-DOES**

[Please contact us for more information and/or a quotation request.](#)

#### Specifications

Brand	Eaton LOFMEM WB
Model	LMWB 30 1.2 DOES
Pore Size $\mu$	1.2

#### Dimensions / Weights

Length in inch	30.0
Length in cm	76.2
Quantity per case	12 pieces



[Download Datasheet](#)

## INFORMACIÓN ADICIONAL

**url\_fab** <https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-30-1.2-DOES/LMWB-30-1.2-DOES-Filter-Cartridge/index.html>

**brand** Lenntech

COTECNO