

## EATON LOFMEM WB LMWB-40-0.2-2S FILTER CARTRIDGE



**SKU:** LMWB-40-0.2-2S | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

COTECNO

## DESCRIPCIÓN DEL PRODUCTO

### EATON LOFMEM WB LMWB-40-0.2-2S FILTER CARTRIDGE

**PART NUMBER: LMWB-40-0.2-2S**

[Please contact us for more information and/or a quotation request.](#)

#### Specifications

Brand	Eaton LOFMEM WB
Model	LMWB 40 0.2 2S
Pore Size $\mu$	0.2

#### Dimensions / Weights

Length in inch	40.0
Length in cm	101.6
Quantity per case	12 pieces



[Download Datasheet](#)

## INFORMACIÓN ADICIONAL

**url\_fab** <https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-40-0.2-2S/LMWB-40-0.2-2S-Filter-Cartridge/index.html>

**brand** Lenntech

COTECNO