

EATON LOFMEM WB LMWB-40-0.3-DOES FILTER CARTRIDGE



SKU: LMWB-40-0.3-DOES | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

DESCRIPCIÓN DEL PRODUCTO

EATON LOFMEM WB LMWB-40-0.3-DOES FILTER CARTRIDGE

PART NUMBER: LMWB-40-0.3-DOES

[Please contact us for more information and/or a quotation request.](#)

Specifications

| | |
|-----------------|------------------|
| Brand | Eaton LOFMEM WB |
| Model | LMWB 40 0.3 DOES |
| Pore Size μ | 0.3 |

Dimensions / Weights

| | |
|-------------------|-----------|
| Length in inch | 40.0 |
| Length in cm | 101.6 |
| Quantity per case | 12 pieces |



[Download Datasheet](#)

INFORMACIÓN ADICIONAL

url_fab <https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-40-0.3-DOES/LMWB-40-0.3-DOES-Filter-Cartridge/index.html>

brand Lenntech

COTECNO