



EATON LOFMEM WB LMWB-40-0.45-2S FILTER CARTRIDGE



SKU: LMWB-40-0.45-2S | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

COTECHNO



DESCRIPCIÓN DEL PRODUCTO

EATON LOFMEM WB LMWB-40-0.45-2S FILTER CARTRIDGE

PART NUMBER: LMWB-40-0.45-2S

[Please contact us for more information and/or a quotation request.](#)

Specifications

Brand	Eaton LOFMEM WB
Model	LMWB 40 0.45 2S
Pore Size μ	0.45

Dimensions / Weights

Length in inch	40.0
Length in cm	101.6
Quantity per case	12 pieces



[Download Datasheet](#)



INFORMACIÓN ADICIONAL

url_fab <https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-40-0.45-2S/LMWB-40-0.45-2S-Filter-Cartridge/index.html>

brand Lenntech

COTECHNO