

EATON LOFMEM WB LMWB-5-0.2-2S FILTER CARTRIDGE



SKU: LMWB-5-0.2-2S | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

DESCRIPCIÓN DEL PRODUCTO

EATON LOFMEM WB LMWB-5-0.2-2S FILTER CARTRIDGE

PART NUMBER: LMWB-5-0.2-2S

[Please contact us for more information and/or a quotation request.](#)

Specifications

Brand	Eaton LOFMEM WB
Model	LMWB 5 0.2 2S
Pore Size μ	0.2

Dimensions / Weights

Length in inch	5.0
Length in cm	12.7
Quantity per case	12 pieces



[Download Datasheet](#)

INFORMACIÓN ADICIONAL

url_fab	https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-5-0.2-2S/LMWB-5-0.2-2S-Filter-Cartridge/index.html
brand	Lenntech

COTECNO