

EATON LOFMEM WB LMWB-5-0.65-3S FILTER CARTRIDGE



SKU: LMWB-5-0.65-3S | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

COTECNO

DESCRIPCIÓN DEL PRODUCTO

EATON LOFMEM WB LMWB-5-0.65-3S FILTER CARTRIDGE

PART NUMBER: LMWB-5-0.65-3S

[Please contact us for more information and/or a quotation request.](#)

Specifications

Brand	Eaton LOFMEM WB
Model	LMWB 5 0.65 3S
Pore Size μ	0.65

Dimensions / Weights

Length in inch	5.0
Length in cm	12.7
Quantity per case	12 pieces



[Download Datasheet](#)

INFORMACIÓN ADICIONAL

url_fab <https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-5-0.65-3S/LMWB-5-0.65-3S-Filter-Cartridge/index.html>

brand Lenntech

COTECNO