

EATON LOFPLEAT AG LPAG-20-0.2-DOES FILTER CARTRIDGE



SKU: LPAG-20-0.2-DOES DuP | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

DESCRIPCIÓN DEL PRODUCTO

EATON LOFPLEAT AG LPAG-20-0.2-DOES FILTER CARTRIDGE

PART NUMBER: LPAG-20-0.2-DOES

[Please contact us for more information and/or a quotation request.](#)

Specifications

Brand	Eaton LOFPLEAT AG
Model	LPAG 20 0.2 DOES
Pore Size μ	0.2

Dimensions / Weights

Length in inch	20.0
Length in cm	50.8
Quantity per case	12 pieces



[Download Datasheet](#)



LOFPLEAT AG SERIES CARTRIDGE FILTERS

Get absolute efficiency with this cost-effective all polypropylene cartridge. The pleated design provides a large surface area for long lasting filtration efficiency.

Features / Benefits

- Fits broad application range with 0.2 to 100 micron retention
- Acceptable for food and beverage applications

- Meets USP Class VI biological tests for plastics
- Beta rating of 5000, 99.98% efficiency
- High flow and long life for minimum maintenance
- Available in lengths up to 40 inches
- Eliminates dirt unloading at high differential pressures due to fixed pore structure

Filter Specifications

Media

- Polypropylene

Inner core

- Polypropylene

End caps

- Polypropylene

Cage

- Polypropylene

Gaskets/O-Rings options

- Buna-N, EPDM, Silicone, Viton®, FEP encapsulated Viton (O-Rings only)

Micron ratings

- 0.2, 0.45, 1.0, 2.5, 5.0, 10, 25, 50, 100 μm

Typical Applications

- process applications
- chemicals

INFORMACIÓN ADICIONAL

url_fab <https://www.lenntech.com/products/Eaton-LOFPLEAT-AG/LPAG-20-0.2-DOES/LPAG-20-0.2-DOES-Filter-Cartridge/index.html>

brand Lenntech

COTECNO