

## EATON LOFPLEAT AG LPAG-20-1-2S FILTER CARTRIDGE



**SKU:** LPAG-20-1-2S DuP | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

## DESCRIPCIÓN DEL PRODUCTO

### EATON LOFPLEAT AG LPAG-20-1-2S FILTER CARTRIDGE

**PART NUMBER: LPAG-20-1-2S**

[Please contact us for more information and/or a quotation request.](#)

#### Specifications

Brand	Eaton LOFPLEAT AG
Model	LPAG 20 1 2S
Pore Size $\mu$	1

#### Dimensions / Weights

Length in inch	20.0
Length in cm	50.8
Quantity per case	12 pieces



[Download Datasheet](#)



## LOFPLEAT AG SERIES CARTRIDGE FILTERS

Get absolute efficiency with this cost-effective all polypropylene cartridge. The pleated design provides a large surface area for long lasting filtration efficiency.

#### Features / Benefits

- Fits broad application range with 0.2 to 100 micron retention
- Acceptable for food and beverage applications

- Meets USP Class VI biological tests for plastics
- Beta rating of 5000, 99.98% efficiency
- High flow and long life for minimum maintenance
- Available in lengths up to 40 inches
- Eliminates dirt unloading at high differential pressures due to fixed pore structure

## **Filter Specifications**

### **Media**

- Polypropylene

#### **Inner core**

- Polypropylene

#### **End caps**

- Polypropylene

#### **Cage**

- Polypropylene

#### **Gaskets/O-Rings options**

- Buna-N, EPDM, Silicone, Viton®, FEP encapsulated Viton (O-Rings only)

#### **Micron ratings**

- 0.2, 0.45, 1.0, 2.5, 5.0, 10, 25, 50, 100  $\mu\text{m}$

#### **Typical Applications**

- process applications
- chemicals

## INFORMACIÓN ADICIONAL

<b>url_fab</b>	<a href="https://www.lenntech.com/products/Eaton-LOFPLEAT-AG/LPAG-20-1-2S/LPAG-20-1-2S-Filter-Cartridge/index.html">https://www.lenntech.com/products/Eaton-LOFPLEAT-AG/LPAG-20-1-2S/LPAG-20-1-2S-Filter-Cartridge/index.html</a>
<b>brand</b>	Lenntech

COTECNO