




## FILTRO CUALITATIVO VELOCIDAD MEDIA 100PCS / CAJA



**SKU:** N / A | **Categorías:** [Equipo Analítico](#) |

## VARIACIONES

Imagen	SKU	Descripción	Capacidad (mm)
 A blue and white rectangular pack of filter paper. The text on the pack reads "FILTER PAPER", "102", and "100 CIRCLES 150mm".	B-01-21-20-000101		110
 A blue and white rectangular pack of filter paper. The text on the pack reads "FILTER PAPER", "102", and "100 CIRCLES 150mm".	B-01-21-20-000102		125
 A blue and white rectangular pack of filter paper. The text on the pack reads "FILTER PAPER", "102", and "100 CIRCLES 150mm".	B-01-21-20-000103		150

## DESCRIPCIÓN DEL PRODUCTO

### Descripción

Filtro cualitativo velocidad media 100pcs / caja

### Especificaciones Técnicas

Código del Producto	Modelo	Capacidad (mm)
B-01-21-20-000101	BellFiltroCual-110	110mm
B-01-21-20-000102	BellFiltroCual-125	125mm
B-01-21-20-000103	BellFiltroCual-150	150mm

## INFORMACIÓN ADICIONAL

COTECNO