


GATOS DE LABORATORIO SS BELL-GATLAB-SS



SKU: B-01-29-0603-0101 | **Categorías:** [Otros equipos de laboratorio](#) |

GALERÍA DE IMÁGENES



DESCRIPCIÓN DEL PRODUCTO

Especificaciones:

| Modelo | Bell-GatLab-SS | Bell-GatLab-AL | Bell-GatLab-R |
|------------|-------------------|-------------------|-------------------|
| Código SKU | B-01-29-0603-0101 | B-01-29-0603-0102 | B-01-29-0603-0103 |
| Talla | 100*100*180mm | 100*100*180mm | 150*150*250mm |
| | 150*150*250mm | 150*150*250mm | 200*200*280mm |
| | 200*200*280mm | 200*200*280mm | |
| | 250*250*330mm | | |
| | 300*300*390mm | | |

INFORMACIÓN ADICIONAL

COTECNO