

MESA VIBRANTE



SKU: B02010126 | **Categorías:** [Sistemas de prueba de mesa vibratoria](#), [Sistemas para pruebas de materiales](#) |

DESCRIPCIÓN DEL PRODUCTO

```
[vc_row type="in_container" full_screen_row_position="middle" column_margin="default" scene_position="center"
text_color="dark" text_align="left" overlay_strength="0.3" shape_divider_position="bottom"
bg_image_animation="none"][vc_column column_padding="no-extra-padding" column_padding_position="all"
background_color_opacity="1" background_hover_color_opacity="1" column_link_target="_self" column_shadow="none"
column_border_radius="none" width="1/1" tablet_width_inherit="default" tablet_text_alignment="default"
phone_text_alignment="default" overlay_strength="0.3" column_border_width="none" column_border_style="solid"
bg_image_animation="none"][vc_gallery type="image_grid" images="129428,129429,129430,129431,129432"
display_title_caption="true" layout="3" item_spacing="4px" gallery_style="1" load_in_animation="none"
img_size="150x150"][tabbed_section style="default" alignment="left" spacing="default" tab_color="Accent-
Color"][/tabbed_section][/vc_column][/vc_row]
```

INFORMACIÓN ADICIONAL

COTECNO