

**PROBETA MEDIDORA DE PLASTICO, BASE HEXAGONAL 250ML, BELLPLA-
PROBINAMPBH-250**



Cod Alt: BellProbPI-250

SKU: B-01-47-0111-0402 B-01-21-0111-0402 | **Categorías:** [Material plástico de laboratorio](#) |

DESCRIPCIÓN DEL PRODUCTO

Modelo	Código Sku	Tamaño disponible ml	Material
BellPla-ProbinaMPBH-50	B-01-47-0111-0400	50	polipropileno
	B-01-21-0111-0400		
BellPla-ProbinaMPBH-100	B-01-47-0111-0401	100	polipropileno
	B-01-21-0111-0401		
BellPla-ProbinaMPBH-250	B-01-47-0111-0402	250	polipropileno
	B-01-21-0111-0402		
BellPla-ProbinaMPBH-500	B-01-47-0111-0403	500	polipropileno
	B-01-21-0111-0403		
BellPla-ProbinaMPBH-1000	B-01-21-0111-0404	1000	polipropileno
BellPla-ProbinaMPBH-2000	B-01-21-0111-0405	2000	polipropileno

INFORMACIÓN ADICIONAL

COTECNO