


## REFRENTADOR DE PROBETAS CILINDRICAS - MARCA COTEST



**SKU:** N / A | **Categorías:** [Equipamiento Metálico COTEST](#) |

## VARIACIONES

Imagen	SKU	Descripción	Medida
 A photograph of a metal test specimen holder for 3-inch cylindrical specimens. It consists of a base plate with a circular opening, a vertical support plate, and a central cylindrical sleeve. The word "cotecno" is printed in small black letters at the bottom right of the image. <small>cotecno</small>	BellRefProC-03 REFRENTADOR DE PROBETAS CILINDRICAS 3"	BellRefProC-03	
 A photograph of a metal test specimen holder for 4-inch cylindrical specimens. It consists of a base plate with a circular opening, a vertical support plate, and a central cylindrical sleeve. The word "cotecno" is printed in small black letters at the bottom right of the image. <small>cotecno</small>	BellRefProC-04 REFRENTADOR DE PROBETAS CILINDRICAS 4".	BellRefProC-04	
 A photograph of a metal test specimen holder for 6-inch cylindrical specimens. It consists of a base plate with a circular opening, a vertical support plate, and a central cylindrical sleeve. The word "cotecno" is printed in small black letters at the bottom right of the image. <small>cotecno</small>	BellRefProC-06 REFRENTADOR DE PROBETAS CILINDRICAS 6".	BellRefProC-06	

## INFORMACIÓN ADICIONAL

COTECNO